VitaVtro Warming Kit

EXPLANATION OF SYMBOLS

The symbol glossary is in line with the ISO 15223-1: Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied-Part 1: General requirements.

Reference number	Symbols	Title
5.1.1		Manufacturer
5.1.2		Authorized representative in the European Community / European Union
5.1.3		Date of manufacture
5.1.4		Use-by date
5.1.5		Batch code
5.1.10		Model number
5.1.11		Country of manufacture
5.2.2		Sterilized using aseptic processing techniques
5.2.8		Do not use if package is damaged and consult instructions for use
5.3.2		Keep away from sunlight
5.3.7		Temperature limit
5.4.2		Do not re-use
5.4.3		Consult instructions for use or consult electronic instructions for use
5.4.6		Contains human blood or plasma derivatives
5.7.7		Medical device
5.7.10		Unique device identifier

EN

INTENDED PURPOSE

VitaVtro Warming Kit is intended for the warming of embryos that have undergone vitrification using VitaVtro Vitrification Kit for ART procedures.

INTENDED PATIENT POPULATION

The devices are applicable to the patient population in need of warming of vitrified embryos.

PRODUCT DESCRIPTION

This product is used for warming vitrified embryos for assisted reproduction technology (ART). It contains three mediums: Human Holding Medium, Human Warming Solution 1, and Human Warming Solution 2.

MODEL NUMBER AND SPECIFICATION

Model No.	HW1	HW2	HHM
JD-1.0	1,5 mL × 2	1,0 mL × 1	1,8 mL × 1
JD-2.0	4,5 mL × 1	4,5 mL × 1	4,5 mL × 2
JD-3.0	4,8 mL × 2	4,8 mL × 1	4,8 mL × 1
JD-4.0	2,0 mL × 5	2,0 mL × 2	2,0 mL × 3
JD-5.0	4,8 mL × 2	4,8 mL × 1	4,8 mL × 2
JD-6.0	2,0 mL × 5	2,0 mL × 3	2,0 mL × 6
JD-7.0	30 mL × 2	30 mL × 1	30 mL × 2

SPECIFICATIONS

Mouse Embryo Assay (1-cell)	≥80% development to blastocyst at 96h
Bacterial Endotoxin	<0.25 EU/mL
pH	7.2-7.6
Sterility	Sterile

LOT-specific test results are available on the Certificate of Analysis (COA) provided with each lot.

COMPONENTS

- Human Holding Medium (HHM) : HEPES Buffered M199 Medium supplemented with Human Serum Albumin (HSA).
- Human Warming Solution 1 (HW1): Same as HHM solution, plus 1.0M sucrose.
- Human Warming Solution 2 (HW2): Same as HHM solution, plus 0.5 M sucrose.
- The reagent vials are made of Polypropylene (PP).

STERILIZATION METHOD

Aseptic processing (including sterilizing filtration and aseptic filling).

STORAGE CONDITIONS

Store at 2°C to 8°C

SHELF LIFE

12 months

PRECAUTIONS

- The user should be a trained professional (e.g. a doctor or embryologist).
- The user should read and understand the user manual and be trained in the correct procedures before using the kit.

WARNINGS

- The long-term safety of embryo cryopreservation is unknown.
- All blood products should be treated as potentially infectious. This product contains Human Serum Albumin (HSA). It was found negative when testing for antibodies to HIV-1/HIV-2, HCV and HTLV-1/HTLV-2 and non-reactive for HBsAg (HBV) and syphilis. However, no known test can guarantee that products derived from humans will not be infectious.
- Not for use in injections.
- Do not use any vials which show evidence of damage, tampering or leaking, or particulate matter and cloudiness.
- To avoid contamination, always use an aseptic technique.
- Do not reuse. After the vial opening, any remaining media should be discarded.
- After use, the product should be disposed of in accordance with relevant regulations of medical waste.

OTHER MATERIALS THAT ARE REQUIRED BUT NOT INCLUDED

- Stereomicroscope
- Heated stage on microscope
- Heated pad close to the microscope
- 1 × four-well dish (4WD)
- 1 × 35 mm Petri Dish
- Stopwatch or timer
- Liquid nitrogen and box
- Pipette

PREPARATION

- All procedures should be performed at room temperature (25-27°C).
- All dishes should be warmed to 25-27°C.
- HHM and HW2 should be warmed to 25-27°C; HW1 should be warmed to 37°C

WARMING AND EXPELLING

The following procedures take the blastocysts for example. The procedures of cleavage embryos are the same.

The volumes of drops in the operation procedures below are recommendatory. Users can adjust the volume according to the results of their laboratory validation.

- Prepare the warming dish (35 mm Petri dish) 1 hour in advance. Pipette 3,000 µL of 25-27 °C HW1 and warm it to 37 °C. The amount of media included in this kit is sufficient to process three blastocysts. A new kit will need to be used if more than this number of blastocysts is to be processed.
- Transfer the liquid nitrogen box containing the storage devices to be warmed to where the warming procedure will be conducted.
- Just before warming, pipette 900 µL of 25-27 °C HW2 solution in well 2 of a Nunc four-well dish. Put 900 µL HHM each into wells 3 and 4.
- Warming of the storage device should be conducted as described in the Instructions for Use for the storage device used.
- The tip of the storage device containing the vitrified blastocysts is immersed in HW1. Ensure that the blastocysts have been recovered from the storage device.

DILUTION AFTER WARMING

- Pipette 20 µL of the HW1 media containing the warmed blastocysts. The blastocysts shall be in the middle of the HW1 (see Figure1). Expel the blastocysts to Well 2 containing HW2. The blastocysts and HW1 should be formed into a sandwich pattern (see Figure 2).


Figure 1. Blastocysts in the middle of HW1


Figure 2. Arrangement of HW1 and HW2 in well 2

- After 3 minutes, pipette 20 µL of the Well 2 solution with the blastocysts to Well 3. Use the sandwich pattern-like step1. Wait 5 minutes. **5 min**
- Transfer the blastocysts to Well 4 with a minimum amount of solution. There is no need to use the sandwich pattern. Wait 5 minutes. **5 min**
- Before performing blastocyst transfer, transfer blastocysts to culture media for 2-4 hours.

Shenzhen VitaVtro Biotech Co., Ltd.
R601, Building B, Hai Ke Xing Tech Park Pingshan District, Baoshan Road No.16 Shenzhen Guangdong, 518118 China
Phone: +86 755 84511813
e-mail: tech@vitavitro.com
EC REP MedNet EC-REP C III GmbH Borkstrasse 10, 48163 Münster, Germany
The electronic IFU (e-IFU) of the device is identical to the hardcopy and can be found together with Summary of Safety and Clinical Performance (SSCP) via the link: www.vitavitro.com
Any serious incident that has occurred in relation to the device could be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.
CE 2797 Doc. No.: IU032207002 Rev. No.: 02A Initial Issue Date: September 5, 2022 Last Revision Date: February 12, 2025

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