



VitaVitro

# Artificial Insemination Catheter



## CAUTION

Federal (U.S.) law restricts this device to sale by or on the order of a physician.

## PRODUCT DESCRIPTION

IUIL-174, IUIS-100. The Artificial Insemination Catheter is composed of a conical joint, an inner tube and an outer tube.

## INDICATION FOR USE

Artificial Insemination Catheter is used for the introduction of washed spermatozoa into the uterine cavity during artificial insemination process.

## CONTRAINDICATION

- 1) Patients with active vaginal or intrauterine infections are prohibited.

- 2) Prohibited for those with sexually transmitted diseases.
- 3) Oviduct factors that cause trouble in sperm and egg bonding.
- 4) Prohibited for pregnant people.

## WARNING

- 1) This product is for the use of assisted reproduction specialists only.
- 2) Operators should be careful during the insemination process to avoid mechanical damage to the vaginal wall and uterine cavity caused by the insemination catheter.
- 3) This product has been sterilized by ethylene oxide, do not use expired.

- 4) This product is sterilized packaging, if damaged, do not use.
- 5) This product is for one-time use only. Please read the instructions carefully before use.
- 6) After use, the product should be disposed of in accordance with relevant regulations of medical waste.

## QUALITY CONTROL TESTING

Sterility tested Passed  
LAL Endotoxin < 5EU/device

## INSTRUCTION FOR USE

- 1) Patients were taken to cross section lithotomy position.
- 2) Exposure of the cervix and uterine opening with a vaginal dilator.

- 3) Clean the cervix with an appropriate cleaning solution.
- 4) Unpack and remove the insemination tube.
- 5) Appropriate amount of sperm washing solution was absorbed with a sterile syringe, connect the syringe to the insemination catheter, push and pull the syringe to wash the insemination catheter, and discard the solution.
- 6) Absorb the appropriate amount of semen, send the insemination catheter slowly through the vagina to the uterine cavity under direct vision, and push the rod of the syringe at a constant speed to inject semen.
- 7) Remove the catheter from the patient.

### DESCRIPTION OF ISO SYMBOLS

The symbol glossary is in line with the SDO-developed standard ANSI/AAMI/ISO 15223-1: Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied-Part 1: General requirements.

Reference number	Symbols	Title of symbol
5.1.1		Manufacturer
5.1.3		Date of manufacture
5.1.4		Use-by date
5.1.5		Batch code

Reference number	Symbols	Title of symbol
5.1.6		Catalog number
5.2.3		Sterilized using ethylene oxide
5.2.8		Do not use if package is damaged
5.3.2		Keep away from sunlight
5.4.2		Do not re-use
5.4.3		Consult instructions for use or consult electronic instructions for use
5.7.7		Medical device



**Shenzhen VitaVibro Biotech Co., Ltd.**

R601, Building B, Hai Ke Xing  
Tech Park, Baoshan Road No.16  
Pingshan District, Shenzhen  
Guangdong, 518118 China

Phone: +86 755 84511813

[www.vitavibro.com](http://www.vitavibro.com)